A propos de la plateforme IRM

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
InroomViewingDevice

For fMRI, patient comfort and instant feedback

Specifications

Screen diagonal	40"
Pixels	3840 x 2160
Contrast	5000:1
Refresh ratio	60Hz @FHD 30Hz@UHD
Camera	1920x1080, 30fps, 1280x720@30fps
Video input	HDMI
Mounting interface	VESA (2 x 100 x 100, D: 400mm, M4)

Multipurpose Usage

The 40" 4K UHD InroomViewingDevice was designed to provide an optimal MR compatible monitor that satisfies the needs of both clinical and advanced scientific applications.

Integrated Camera

With its slim design, high definition display and superior image quality, the InroomViewingDevice is an optimal choice for an easy to use alternative to conventional projectors or goggle based image delivery systems. The innovative, front-facing camera provides an uninterrupted patient surveillance during examination. Thanks to the built-in USB hub, connecting patient communication and interface device is now extremely easy.

Flexible Positioning

The low weight and height adjustable mobile foot stand allows easy positioning of the monitor anywhere in the MRI room. The monitor can also be ceiling mounted (through 3rd party).

Instant Feedback

The monitor facilitates the examination process by allowing the operating personnel to remain inside the examination room during procedures, thus allowing uninterrupted patient care and quick response time, which significantly improves clinical workflow.

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**Pyka 4 Button Handheld
Response Grip**

**(Cambridge research)
(NordicNeuroLab)**

2 Button Bimanual

(Current Design)

.IRM préclinique 7T de MR Solutions

- Antennes : Mouse head coil, Mouse Body coil, Marmoset head coil, Rat head coil
- Logiciels : Precinical Scan, Power scan.
- Matériel de monitoring et d'anesthésie : Minerve Veterinaire

.Simulateur IRM

- MoTrak® Head Motion Tracking System

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MAGNETIC RESONANCE (MR) PROCEDURE **SCREENING FORM FOR VOLUNTEERS**

Last name : _____

First name : _____

Date of Birth : _____
dd/mm/yyyy

Sex : Male () Female () Other ()

Height : _____

Weight : _____

Have you had prior surgery or an operation ?

	Yes	No	If yes, indicate the type and date
Head	<input type="checkbox"/>	<input type="checkbox"/>	
Heart	<input type="checkbox"/>	<input type="checkbox"/>	
Abdomen or Chest	<input type="checkbox"/>	<input type="checkbox"/>	
Orthopaedic surgery	<input type="checkbox"/>	<input type="checkbox"/>	
Vertebral column - spine	<input type="checkbox"/>	<input type="checkbox"/>	
Other	<input type="checkbox"/>	<input type="checkbox"/>	

Please indicate if you have any of the following :

	Yes	No	
Aneurysm clip(s)	<input type="checkbox"/>	<input type="checkbox"/>	
Cardiac pacemaker	<input type="checkbox"/>	<input type="checkbox"/>	
Heart valve prosthesis	<input type="checkbox"/>	<input type="checkbox"/>	
Implanted cardioverter defibrillator (ICD)	<input type="checkbox"/>	<input type="checkbox"/>	
Metallic stent, filter, or coil	<input type="checkbox"/>	<input type="checkbox"/>	
Electronic or Magnetically-activated implant or device	<input type="checkbox"/>	<input type="checkbox"/>	
Neurostimulation system	<input type="checkbox"/>	<input type="checkbox"/>	
Spinal cord stimulator	<input type="checkbox"/>	<input type="checkbox"/>	
Internal electrodes or wires	<input type="checkbox"/>	<input type="checkbox"/>	
Cochlear, otologic, or other ear implant	<input type="checkbox"/>	<input type="checkbox"/>	
Insulin or other infusion pump or device	<input type="checkbox"/>	<input type="checkbox"/>	
Any type of prosthesis (eye, penile, etc.)	<input type="checkbox"/>	<input type="checkbox"/>	
Shunt (spinal or intraventricular)	<input type="checkbox"/>	<input type="checkbox"/>	
Vascular access port and/or catheter	<input type="checkbox"/>	<input type="checkbox"/>	
Medication patch (Nicotine, Nitroglycerine)	<input type="checkbox"/>	<input type="checkbox"/>	
Any metallic fragment or foreign body	<input type="checkbox"/>	<input type="checkbox"/>	
Tissue expander (e.g., breast)	<input type="checkbox"/>	<input type="checkbox"/>	
Surgical staples, clips, or metallic sutures	<input type="checkbox"/>	<input type="checkbox"/>	
IUD, diaphragm, or pessary	<input type="checkbox"/>	<input type="checkbox"/>	
Tattoo or permanent makeup	<input type="checkbox"/>	<input type="checkbox"/>	
Body piercing jewellery	<input type="checkbox"/>	<input type="checkbox"/>	
Hearing aid	<input type="checkbox"/>	<input type="checkbox"/>	
Breathing problem or motion disorder	<input type="checkbox"/>	<input type="checkbox"/>	
Other implant	<input type="checkbox"/>	<input type="checkbox"/>	

	Yes	No
Have you had a prior MRI examination or MR procedure ?	<input type="checkbox"/>	<input type="checkbox"/>
Have you ever been injured by a metallic object or foreign body ?	<input type="checkbox"/>	<input type="checkbox"/>
For female volunteers		
Are you pregnant or experiencing a late menstrual period?	<input type="checkbox"/>	<input type="checkbox"/>
Are you taking any type of fertility medication or having fertility treatments?	<input type="checkbox"/>	<input type="checkbox"/>

IMPORTANT INSTRUCTIONS

- Any operation performed less than 6 weeks before the scheduled date of the MRI exam is a contraindication to the performance of the MRI exam.
- Depending on the operation, even if beyond 6 weeks, the operating protocol may be required.
- Before entering the MR environment or MR system room, you must remove all metallic objects including hearing aids, dentures, partial plates, keys, beeper, cell phone, eyeglasses, hair pins, barrettes, jewellery, body piercing jewellery, watch, safety pins, paperclips, money clip, credit cards, bank cards, magnetic strip cards, coins, pens, pocket knife, nail clipper, tools, clothing with metal fasteners, & clothing with metallic threads.
- Please consult the MRI Technologist or Research staff if you have any question or concern BEFORE you enter the MR system room.

I attest that the above information is correct to the best of my knowledge. I read and understood the contents of this form and had the opportunity to ask questions regarding the information on this form and regarding the MR procedure that I am about to undergo.

Signature of Person Completing Form : _____ **Date :** _____

Form Information Reviewed By (Research staff) :

Name & Signature _____ **Date :** _____

Form Information Reviewed By (MRI Technologist):

Name & Signature _____ **Date :** _____